Viv Highback Chair Technical Specifications

Testing

Our chairs are designed to meet or exceed all EN performance requirements per **EN16139:2013**, which includes the following tests:

- Seat and back static load test
- Seat front edge static load test
- Vertical static load on back
- Foot rest and leg rest static load test
- Arm sideways static load test
- Arm downwards static load test
- Vertical upwards static load on arm rests

- Seat and back durability test
- Seat front edge durability test
- Leg forward static load test
- Leg sideways static load test
- Seat impact test
- Back impact test
- Arm impact test

Structural Testing Certifications

VIVHCH5SL

Viv Highback Chair 5 star with lift mechanism **BS EN 16139,2013, Level 2** Static load 200kg

VIVHACH5SL

Viv Highback Armchair on 5 star with lift mechanism BS EN 16139,2013, Level 1 Static Load 160kg ANSI/BIFMA X5.4 - 2012 Static Load 300lbs

VIVHLO4S

Viv Highback Lounge on 4 star base BS EN 16139,2013, Level 2 Static Load 200kg

Flammability

The cold cure foam material has been tested to meet the following standards:

British standard 5852:2006 Part 2

California Technical Bulletin **117-2013**

Italian standards
UNI 9175-CLASS 1.IM

The flammability standard of the covering fabric is dependent on the choice made by the customer, details of the flammability rating of specific fabrics can be found on supplier's websites listed below.

www.kvadrat.dk
www.camirafabrics.com
www.butefabrics.com
www.panaz.com
www.maharam.com
www.vescom.com
www.ultrafabricsllc.com
www.crestleather.com

Servicing



If a component is damaged and needs replacing then please email our team with an image of the damaged part and reference your order number, we will be happy to facilitate any replacements.

Warranty



naughtone offers a 10 Year warranty across all products Details on the terms and conditions of naughtones warranty can be found at naughtone.com/warranty

Viv Highback Chair Technical Specifications











and I		uct Code Description ensions quoted are rounded to rest 5mm or nearest 1/2 inch	A Overall Width	B Overall Depth	C Overall Height	D Seat Height	E Seat Width	F Arm height	G Height from Seat to Top	H Weight
	1	VIVHLO4S Viv Highback Lounge Chair 4 Star Base	620 mm 24.5 in	640 mm 25 in	930 mm 36.5 in	400 mm 15.5 in	525 mm 20.5 in	N/A	575 mm 22.5 in	18.3kg 40.4lbs
	2	VIVHALO4S Viv Highback Lounge Armchair 4 Star Base	620 mm 24.5 in	640 mm 25 in	930 mm 36.5 in	400 mm 15.5 in	525 mm 20.5 in	175 mm 7 in	575 mm 22.5 in	20.2kg 44.5kg
	3	VIVHCH4SL Viv Highback Chair 4 Star Base with Lift Mechanism	620 mm 24.5 in	640 mm 25 in	1025-1125 mm 40-44 in	450-550 mm 17.5-21.5 in	525 mm 20.5 in	N/A	575 mm 22.5 in	17.5kg 38.6lbs
	4	VIVHACH4SL Viv Highback Armchair 4 Star Base with Lift Mechanism	620 mm 24.5 in	640 mm 25 in	1025-1125 mm 40-44 in	450 – 550 mm 17.5 – 21.5 in	525 mm 20.5 in	175 mm 7 in	575 mm 22.5 in	19.3kg 42.5lbs
	5	VIVHCH5SL Viv Highback Chair 5 Star Base with Lift Mechanism	635 mm 25 in	640 mm 25 in	1025-1125 mm 40-44 in	450-550 mm 17.5-21.5 in	525 mm 20.5 in	N/A	575 mm 22.5 in	18.0kg 39.7lbs
	6	VIVHACH5SL Viv Highback Armchair 5 Star Base with Lift Mechanism	635 mm 25 in	640 mm 25 in	1025-1125 mm 40-44 in	450 – 550 mm 17.5 – 21.5 in	525 mm 20.5 in	175 mm 7 in	575 mm 22.5 in	19.8kg 43.7lbs

